


SARS-CoV-2 Antigen Test

(Lateral Flow Method)



- Easy to use, no equipment required , no PCR Lab required
- Instant result in 15 minutes
- Room temperature storage (2~30°C)
- Noninvasive sampling for Nasopharyngeal or oropharyngeal swab
- Direct detection of virus , faster than PCR test

Operation procedure



Order information

Catalog No.	Product Name	Sample Type	Storage Condition	Shelf Life	Qualification
W196	SARS-CoV-2 Antigen Test (Lateral Flow Method)	Nasopharyngeal swab or oropharyngeal swab	2~30 °C	12 Months	CE

